

SBS-offline for SBSGenericDetector and BBShower

July 12, 2021

1 Introduction

1.1 SBSGenericDetector

SBSGenericDetector is a THaNonTrackingDetector. For the SBSGenericDetector, a detector can be classified as using ADC and/or TDC information. The ADC can be different modes: kNone (no ADC), kADCSimple (like a 792 or fastbus), kADC (F250 in mode 7) and kWaveform (F250 in mode 1). ADC data from a F250 in mode 7 can have up three pulses recorded per channel per event. ADC data from a F250 in mode 1 is a vector of 4ns samples. The size of the vector is fixed in the configuration of the FADC readout list. The size of the vector is determined by the FADC250 module decoder. The TDC can be in modes: kNone (no TDC for detector), kTDCSimple (only leading edge), kTDC (leading and trailing edge). The constructor sets variables to default values which are given in Table 1. Some variables can be set by functions that can be called from your replay script. You can also set this in the constructor of the detector that inherits from SBSGenericDetector. If one sets SBSModeADC::Mode not equal to kADC then SetDisableRefADC set to false.

1.2 SBSCalorimeter and SBSBBShower

SBSCalorimeter constructor inherits from the SBSGenericDetector and leaves the SBSModeADC::Mode as kADCSimple and SBSModeTDC::Mode as kNone. SBSBBShower constructor inherits from the SBSCalorimeter and sets SetModeADC(SBSModeADC::kWaveform) and SetModeTDC(SBSModeTDC::kNone).

Variable	Default	Comment
fNrows	0	
fNcolsMax	0	
fNlayers	0	
fModeADC	kADCSimple	SetModeADC(SBSModeADC::Mode mode)
fModeTDC	kNone	SetModeTDC(SBSModeTDC::Mode mode)
fDisableRefADC	true	SetDisableRefADC(Bool_t b)
fDisableRefTDC	true	SetDisableRefTDC(Bool_t b)
fConst	0	
fSlope	0	
fAccCharge	0	
fStoreRawHits	false	SetStoreRawHits(Bool_t var)
fStoreEmptyElements	true	SetStoreEmptyElements(Bool_t b)

Table 1: Table of variables in the constructor

2 SBSGenericDetector ReadDatabase

A detector is defined as a grouping of elements which are arranged in rows, cols and layers. The element is defined by the SBSElement class which stores information about the element such as its geometry and the TDC and/or ADC data. The TDC and ADC data is stored in structures and classes defined in the SBSData class. There are two numbering schemes for vectors of the elements. Either by element number and like a 3-d array with row, col and layer. ncols is a parameter that needs to be set. nrows is optional and the default is 1. ncols is an array. If the size of ncols array is 1, then all the nrow are given the same number of columns equal to ncols[0]. If the ncols array is nrow size, then each row has a different number of columns. The nlayers is optional and the default is 1.

Need to set the vector parameters xyz (the x,y and z starting position for the geometry) and dxyz (the dx, dy and dz spacing between elements).

The parameter start_chanmap is optional and is not used in the code. The parameter row_offset_pattern is optional and is not used in the code. Maybe these can be eliminated.

The detmap gives list of the crate, slot, start channel, end channel and reference time index for each module that is used by the detector. One can also add the module number in the detmap, if one sets the parameter model_in_detmap = 1. By default model_in_detmap = 0 and setting the de-

tector map is probably only necessary for when using the VETROC. When one wants to include the model in the detector map, then the detmap format changes and should be crate, slot, start channel , end channel, model number and reference time index for each module that is used by the detector. The reference time index is an integer that indicates which reference time indicated in the chanmap to use for that module. The detmap vector is an input to the FillDetMap function of the THaDetMap which produces the object fDetMap.

The chanmap is vector of the element index for the crate, slot and channel. The chanmap order must follow the ordering of the crate and slots in the detmap. It is an optional parameter. Without the chanmap defined, the default is to assign the elements to each channel in the detmap in order. When creating the chanmap, the element index starts at zero and the maximum is the total number of elements (nrow x ncol x nlayer) minus one. The chanmap format is to match the number of channels in the detmap and in the same order of the detmap. If a channel is not used by the detector then put a -1 in the chanmap. If the channel is a reference time, then put a -1000 in the chanmap.

The code loops through the modules in the fDetMap to create a fChanMap and fRefChanMap based on the chanmap. fChanMap and fRefChanMap are use in Decode method to determine which channel info is loaded into which element. This determines the number of elements, fNelem, and the number of reference time elements for the TDC, NRefADCElem and ADC, NRefADCElem. The vector size of fRefChanMap is equal to the number of -1000 channels that are found and the index of the fRefChanMap vector is the order that it finds the -1000. So in the detmap , the reference time index refers to the order of the -1000 in the chanmap and goes from 0 to the number of chanmap =-1000 minus one. As an example for the FADC, there can be one reference time in one module (reference time index = 0) which is used by all FADC modules (so detmap would have reference time index = 0 for each FADC). Similarly for the TDC, there can be one reference time in one module to be used for all modules or a reference time for each module. For the case of FADC and TDC modules, if the TDC modules are first in the detmap and have one module with chanmap = -1000 then the TDC detmap would use reference time index = 0 and if the FADC modules with chanmap = -1000 then the ADC detmap would use reference time index = 1.

The code then reads in parameters used to process the TDC or ADC raw data into calibrated data for each element. All these parameters are optional with default values in the code. For the TDC, the parameters are given in Table 2. The parameter vector can have all elements set individually or

have one value which is used for all elements. The same parameters can be read-in for the TDC reference time but the "tdc" is replaced by "reftdc".

Parameter name	Description	Global variable
tdc.offset	offset (default=0)	tdc_offset (unit channels)
tdc.cal	calibration (default=0)	tdc_cal (ns/channel)
tdc.GoodTimeCut	Expected Good Time peak (ns) (default=0)	tdc_GoodTimeCut

Table 2: Table of TDC parameters

For the ADC, the parameters for the kADCSimple are given in Table 3. For processing the kWaveform ADC data to extract the pulse information ,

Parameter name	Description	Global variable
adc.pedestal	Pedestals (not used for waveform)	adc_ped
adc.gain	gain to convert to GeV	adc_gain

Table 3: Table of kADCSimple parameters

there are the following parameters given in Table 4 are used in addition to the kADCSimple parameters. The parameter vectors can have all elements set individually or have one value which is used for all elements. The same parameters can be read-in for the ADC reference time but the "adc" is replaced by "refadc".

adc.conv	conversion factors (waveform is mV/chan)	adc_conv
adc.thres	threshold (mV)	adc_thres
adc.FixThresBin	fixed threshold crossing bin	adc_FixThresBin
adc.NSB	# of integration samples before threshold crossing	adc_NSB
adc.NSA	# of integration samples after threshold crossing	adc_NSA
adc.NPedBin	# of samples for pedestal determination	adc_NPedBin

Table 4: Table of ADC parameters

The code creates a vector of elements called fElements based on the SBSElement class. The element number is incremented by looping through number of columns inside the looping through rows. The element is created and loads the x,y,z locations along with the row, column, layer and element number. The element can be have ADC and/or TDC information depending

on how the GenericDetector is defined. If the GenericDetector has TDC info, then the element is set to have an SBSData::TDC data function and the offset, calibration and GoodTimeCut are loaded. If the GenericDetector has ADC info, then it depends on the ADC mode. If the ADC mode is kWaveform, the element is set to have an SBSData::Waveform data function and the pedestal (which is not used), gain and conversion factor are loaded. For the kWaveform there is an additional call to load the parameter used to determine ADC signal in the waveform. If the ADC mode is not kWaveform, the element is set to have an SBSData::ADC and the pedestal and gain are loaded.

The code creates a vector of reference time elements called fRefElements bases on the SBSElement class. It follows the same logic for creating the reference time elements as the detector elements, except that the reference time element will either be a TDC or an ADC reference time element. When the reference time element is created the row number and element number are both the reference time index.

2.1 BBShower

The BB Preshower is 26 rows of 2 columns in one layer. The BB SHower is 27 rows of 7 columns in one layer. For the BB Preshower, the left side was put in slot 3 and slot 4 (chan 9) and the right side slot 3 (chan 10) to slot 6 (chan 3). The left side is column 0, so all the even element indices. The right side is column 1, so all the odd element indices. For BB Shower, the signals are in slot 6 (chan 4) to slot 20 (chan 0). As an example, row 0 and col 0-6 is in slot 6 and channels 4 to 10.

3 SBSGenericDetector Decode

First the ClearEvent method is called. This clears the ADC and/or TDC data in all the elements.

If the reference times are in the channel map, then the code looks for the reference time in the data. The code loops the channels in each module. If there is event data for the reference time channel, then the data is processed for the FRefElement with DecodeADC or DecodeTDC depending if the module is an ADC or TDC. Only expect ADC reference times for the FADCs. For ADC waveform, the time for the pulse in the reference time is determined the same way as the time a regular detector element. The TDC reference times are processed the same as a regular detector element.

The total number of ADC and TDC reference channels hit in the Decode method is fNRefhits.

The code loops through the channels in each module. If there is event data for the detector channel, then the data is processed for that fElement with DecodeADC or DecodeTDC depending if the module is an ADC or TDC. The total number of ADC and TDC channels hit in the Decode method is fNhits. The ADC and TDC channels are not counted separately. **Look into adding counter for the ADC and TDC separately and put in the tree.**

3.1 SBSGenericDetector DecodeTDC

The method assumes that the reference time elements will be processed first in the Decode method.

If a detector time is being processed, then the method first determines the good reference time to use for the detector element based on the module's reference time index. It loops through the reference time TDC hits and selects the reference time hit that has the smallest $\text{abs}(\text{Tdc} - \text{GoodTimeCut})$. The GoodTimeCut is in ns. This good reference time is the raw time in units of channels and is locally stored in the variable reftime. If a reference time is being processed, then reftime is set to zero.

For either a reference or detector time, the method loops through the hits for that element and first determines whether the hit is leading (LE) or trailing edge (TE). If the first hit is a TE, then it skips that hit. If the last hit is a LE, then that it skips that hit. If it has passed the first two tests, then the hit is processed using the SBSData method TDC::Process. TDC::Process has the arguments of the element ID, raw time value and edge value (LE = 0, TE = 1). The raw time value that is loaded is the raw time value from the channel minus reftime (local variable determined earlier in DecodeTDC).

TDC::Process fills the fTDC which is a SBSData::TDCData container for that detector or reference time element. TDCData holds the offset, calibration, GoodTimeCut and the good_hit. The offset, calibration, GoodTimeCut are set in SBSGenericDetector::ReadDatabase. good_hit is an integer variable which holds the index of the good TDC hit that is determined by the SBSGenericDetector::FindGoodHit method that is described later. TDCData also contains a vector, "hits" of the SBSData::TDCHit structure. The SBSData::TDCHit structure contains the "le" SingleData structure, "te" SingleData structure, "ToT" SingleData structure and "elem ID", the integer of the element ID. The SBSData::SingleData structure contains the

"raw" and "cal" (the raw and calibrated values). For a TDC, the raw value is reference time subtracted and in units of channels and the calibrated value is the (raw value - offset)*cal.

TDC::Process has a counter array fEdgeIdx[edge] with edge = 0 (LE) or 1 (TE). If the TDC detects LE and TE, then the goal is to first find a LE hit and match it with a TE hit which usually follows the LE. If the counter fEdgeIdx[edge] is larger or equal to the size of the TDCHit hits vector for that element, then TDC::Process increments the size of the hits vector and the fEdgeIdx[edge] counter is incremented. The element ID is loaded into the hits vector. If the edge is LE then the "le" structure is filled. If the edge is TE then the "te" structure is filled. If the fEdgeIdx[0] equals fEdgeIdx[1] then the "ToT" structure is filled with the "te" - "le". If the fEdgeIdx[1] is larger than fEdgeIdx[0], then fEdgeIdx[0] is set equal to fEdgeIdx[1]. This is done to insure that if the next hit is a LE then the size of the hits vector will be increased whether the next hit is a LE or TE. The TDCHit hits vector is ideally assuming that it will have an LE and TE (or just LE if the TDC module was programmed that way). If there is only a LE or TE filled by a hit, then the missing information will have a value of zero. **Should the default value be kBig?**

After looping through all the hits for that element, the element will have an TDC data structure that contains the TDCHit hits vector. In the SBSGenericDetector::CoarseProcess, there is a call to SBSGenericDetector::FindGoodHit for each element that has TDC data and the hit which is the "good" hit is chosen.

3.2 SBSGenericDetector DecodeADC

If ADC mode is kADCSimple, the SBSData::ADC Process method is called with a single argument. It fills the fADC which is a SBSData::ADCData container for that detector or reference time element. The SBSData ADCData structure holds the pedestal, calibration, time calibration (for FADC) and the amplitude calibration. In addition, ADCData contains vector of the PulseADCData structure, hits, and the good_hit which is the index number of the good hit which are needed for multihit FADC. For a kADCSimple mode , only the vector will just be the single hit. The PulseADCData is a structure containing the integral, time and amplitude which are all SingleData structures having the "raw" and "val" (calibrated value). For the kADCSimple mode, the ADC::Process with single argument fills the SingleData integral with the raw data and the (raw-ped)*cal with the time and amplitude set to zero and sets fHasData to true.

The ADC mode kADC which is for mode = 7 in the FADC250 needs to be worked on

When ADC is in waveform mode, a sample vector is filled with ADC value for each 4ns bin. The SBSData::Waveform::Process is called. A vector with the raw ADC converted to mV is filled. The pedestal (in mV) is determined by averaging the first NPedBin ADC values. A vector with the pedestal subtracted ADC is filled. The threshold crossing bin is found by looping through the sample raw ADC vector and finding the first bin which has a raw ADC value larger than the pedestal plus the threshold. If a threshold crossing bin, TC, is found, then the ADC samples are integrated from TC-NSB to TC+NSA-1. If no threshold crossing bin, TC, is found, then the ADC samples are integrated from FixThresBin-NSB to FixThresBin+NSA-1. The ADC integration is converted to pC using the 4ns bin width and 50 ohm resistance.

If a threshold crossing bin is found, then the peak amplitude is found within the integration region. The peak, Vpeak, is the ADC value in the sample which within the integration window (starting from the threshold crossing) the following bin has a smaller ADC value than previous sample. To calculate the pulse time, first $V - Mid = (V_{peak} + pedestal)/2$ is calculated. Then loops through the integration region and finds the sample, i, with which has value greater than Vmid and the next sample, i+1, is less than Vmid. The time is $4ns * [(i + (V_{Mid} - V_i)/(V_{i+1} - V_i))]$. If no threshold crossing bin is found, then the time and amplitude are set to zero. For now the code only find the first pulse that is above threshold.

3.3 SBSGenericDetector::CoarseProcess

The CoarseProcess method separately loops through the reference time elements and the detector elements. The goal is to fill the fRefGood and the fGood SBSGenericOutputData structures. The SBSGenericOutputData is a structure the contains vectors of the row,col,layer, element ID. SBSGenericOutputData contains vectors of the ADC integral, pedestal subtracted ADC integral, the calibrated pedestal subtracted ADC integral, FADC pulse amplitude, FADC pedestal subtracted pulse amplitude and the FADC pulse time. SBSGenericOutputData contains vectors of the TDC LE time, TE time and Tot (TE-LE) time. SBSGenericOutputData contains vectors of number of samples, index for the start of samples and the samples of the full FADC waveform.

The method first loops through the reference time elements. The fRefGood structure is filled with row, column, layer and element ID. If the refer-

ence time element is a TDC and has TDC data, then the TDCHit structure for the good reference time hit is retrieved and the fRefGood TDC structure is filled. By default SBSGenericDetector sets fStoreEmptyElements to true, but the detector class that inherits could set fStoreEmptyElements to false or this can be set in the replay script. If fStoreEmptyElements is true then the TDC structure is filled with kBig.

If the element has TDC data then it is possible have the code store the TDC hit information for all hits in the reference time elements in the fRefRaw SBSGenericOutputData structure. By default SBSGenericDetector sets fStoreRawHits to false but the detector class that inherits could set fStoreRawHits to true or this can be set in the replay script. If true, then the fRefRaw structure is filled with the element ID and the calibrated LE, TE and Tot values for all hits in an element for all elements in the loop that have TDC data.

Need to add code for processing the ADC reference times

The method then loops through all the detector elements. **It calls SBSElement::CoarseProcess which fills the energy for the element. This should be eliminated?.** The SBSGenericDetector::FindGoodHit method is called for the element. If the element has TDC data then the "good" TDC hit is determined by finding the TDC hit which has the smallest abs(Tdc - GoodTimeCut). The GoodTimeCut is in ns. **If the element has ADC data, then first hit is used. Presently no multihit with the waveform** The detector class that inherit from GenericDetector could implement their own FindGoodHit method.

After the good hit is determined, if the element has TDC data, then fGood is filled with TDC information. If fStoreEmptyElements is true then the TDC information is filled with kBig. If fStoreRawHits is true, then all the hits in the element is loaded into the fRaw SBSGenericOutputData structure.

Next, if the element has ADC data (fHasData is true), then the "good" ADC hit is used to fill the fGood with ADC information. If the ADCmode is kADCsimple, then the fGood pedestal, raw integral ("a"), raw -pedestal ("a_p") and the calibrated value ("a_c") are loaded. For kADC mode, the amplitude ("a_amp"), amplitude - pedestal ("a_amp_p") and pulse time ("a_time") are loaded. **need to work on this but at this point noone uses this adc mode.** If fStoreRawHits is true (by default it is false) then the calibrated integral, calibrated amplitude and pulse time are loaded into fRaw structure. **Need to add the ADC element ID, since it will be different than the TDC element ID** If fStoreEmptyElements is true (it is set true by default) then the ADC information is filled with zero. If the

ADC mode is Waveform, then the fGood pedestal, raw integral ("a"), raw -pedestal ("a_p"), the calibrated value ("a_c"), the amplitude ("a_amp"), amplitude - pedestal ("a_amp_p") and pulse time ("a_time") are loaded. In addition, the index of the sample data in the "samps" vector ("sidx"), the number of samples per element ("nsamps") and the calibrated sample data ("samps") are loaded into a vector for all elements. **Presently, there is no check on a "good hit", since there is only one hit determined per waveform.**

In general, a detector class that inherits from the GenericDetector should first call GenericDetector::CoarseProcess and then use the fGood structure to process the data for further analysis.

3.4 SBSGenericDetector Tree variables

For the tree variables, most are arrays of variable size depending on the number of hits. The total number of ADC and TDC channels hit in the Decode method is fNHits. The ADC and TDC channels are not counted separately.

Tree name	Description	Global variable
nhits	Number of total hits	fNHits
ngoodhits	Number of total good hits	fNGoodhits

If fStorerawHits is set true, then the ADC/TDC raw info is stored.

Tree name	Description	Global variable
hits.t	all hit's Calibrated Leading edge TDC	fRaw.t
hits.t_te	all hit's Calibrated Trailing edge TDC	fRaw.t_te
hits.t_tot	all hit's TDC Time Over Threshold	fRaw.t_ToT

Tree name	Description	Global variable
hits.a	all ADC integral of the all hits (units pC)	fRaw.a
hits.a_amp	all ADC peaks of the all hits (units mV)	fRaw.a_amp
hits.a_time	all ADC time of the all hits (units mV)	fRaw.a_time

The good hit information is the following variables.

Tree name	Description	Global variable
nhits	Number of total hits	fNHits
ngoodhits	Number of total good hits	fNGoodhits
row	Row of the detector element	fGood.row
col	Column of the detector element	fGood.col
layer	Layer of the detector element	fGood.layer
ped	Pedestal of the good hit (units mV)	fGood.ped

Tree name	Description	Global variable
a	ADC integral of the good hit (units pC)	fGood.a
a_p	Pedestal subtracted ADC integral of the good hit (units pC)	fGood.a_p
a_c	Energy Calibrated Pedestal subtracted ADC integral of the good hit (units GeV)	fGood.a_c

Tree name	Description	Global variable
a_amp	ADC peak of the good hit(units mV)	fGood.a_amp
a_amp_p	Pedestal subtracted ADC peak of the good hit (units mV)	fGood.a_amp_p
a_time	ADC time of the good hit (units ns)	fGood.a_time
a_amp	ADC peak of the good hit (units mV)	fGood.a_amp

Tree name	Description	Global variable
tdc	Calibrated Good hit Leading edge TDC	fGood.t
tdc_te	Calibrated Good hit Trailing Edge TDC	fGood.t_te
tdc_tot	Calibrated Good hit Time over Threshold	fGood.t_Tot

Tree name	Description	Global variable
samp_idx	Index in the samples vector for given row-col element	fGood.sidx
samp_nsamps	Number of samples for a given row-col element	fGood.nsamps
samp	Calibrated Pedestal subtracted ADC samples (units mV) in 4ns bins	fGood.samps

4 Short description of FADC250

The FADC250 is a fast data acquisition electronics module which is able to sample a signal pulse at 250MHz (i.e. 4ns time sample). It has two main data acquisition modes:

- Mode 1: all samples within a defined time window are recorded; this allows to perform the waveform reconstruction offline.
- Mode 7: if the recorded pulse passes a defined threshold, the following quantities are recorded: the pedestal, the pulse time, the pulse amplitude, and the pulse integral over a defined time window; Up to 3 pulses can be recorded the trigger window

For the SBS hadron calorimeter (HCal) and SBS BBSHower/PreSHower, the FADC will be used with acquisition mode 1.

Details on the algorithm of the FADC is given in the manual. Here is a short description. In mode 7, a threshold is defined in the ROL FADC configuration file for each channel in the FADC. When a pulse has one sample that passes the threshold, then that sample is marked as the threshold crossing sample, TC. The pulse is integrated from NSB samples before the TC and NSA samples after the TC. NSB and NSA are programmed in the ROL FADC configuration file. The average of first four samples in the trigger window are used to determine the pedestal or VMIN. The peak amplitude (VPEAK) is determined by finding a sample beyond TC for which the sample value first decreases. The fine time, TF, is determined in bins of 62.5ps. First the time bin is roughly determined by the time bin of the half amplitude ($VMID = (VPEAK + VMIN) / 2$). Then the fine time is determined as $TF = 64 * (VMID - V(N1)) / (V(N1+1) - V(N1))$. The sample number N1 is found on the leading edge of the pulse that satisfies $V(N1) \leq VMID < V(N1+1)$.